



SMOBIO

Small Bio, Smart Tool

www.smobio.com

Product Information

GetClone™ PCR Cloning Vector

CV1000	20 RXN		
	pGet1.1 Vector (50 ng/μl)	23 μl	
	pGet-For Primer (10 μM)	100 μl	
	pGet-Rev Primer (10 μM)	100 μl	

Storage

-20°C for 24 months

Features

1. Cloning efficiency greater than 90%
2. IPTG and X-Gal are not required
3. Accepts a wide range of insert/vector ratios 0.5:1 to 12:1
4. Accepts insert size from 6 bp to 11 kb
5. The phosphorylation of PCR fragments is not required
6. Accepts **blunt end** amplicon or DNA fragment (not for sticky ends)

Description

The GetClone™ PCR Cloning Vector is a positive selection system for high efficiency cloning of blunt end DNA or PCR products amplified by high fidelity DNA polymerase. The linearized GetClone™ pGet1.1 Vector contains a lethal gene which can be disrupted by ligation of a blunt end DNA insert into the cloning site. After ligation and transformation, only *E.coli* clones carrying the pGet1.1 Vector with inserted DNA at cloning site are able to propagate on LB-ampicillin agar plates, eliminating the additional needs of IPTG and X-Gal for blue/white screening.

Primers Sequence

pGet-For Primer:

5'-TCGAAGTTAAAGATGATTACGG-3'

pGet-Rev Primer:

5'-TCTCTCGATAGCATTTCCTGC-3'

Ligation Example 1 (NEB T4 DNA Ligase #M0202)

Insert (Blunt end)	X μ l (Y ng*)
pGet1.1 (2995 bp)	1 μ l (50 ng)

Mix well then add

10X T4 DNA Ligase Buffer	2 μ l
T4 DNA Ligase	1 μ l
ddH ₂ O	to 20 μ l
Final volume	20 μ l

Mix well then incubate at 16°C or room temperature (20~25°C) for 1 hours.

Ligation Example 2 (TOYOBO Ligation High ver2 #LGK-201)

Insert (Blunt end)	X μ l (Y ng*)
pGet1.1 (2995 bp)	1 μ l (50 ng)
ddH ₂ O	up to 7 μ l
Ligation high ver2	3.5 μ l
Final volume	10.5 μ l

Mix well then incubate at 16°C or room temperature (20~25°C) for 5~30 mins.

*For 3/1 of Insert/Vector molar ratio:

$$Y(\text{ng}) = \frac{3}{1} \times \frac{50(\text{ng}) \times \text{Insert size (kb)}}{2.995 (\text{kb})}$$

Transformation

The GetClone™ is compatible with most available competent *E. coli* cells. Apply 1 ~10 µl of ligation mixture to 10 times volume competent *E. coli* cells. Perform transformation procedures according to the instruction of the competent cells. Spread the transformed *E. coli* cells on an LB-ampicillin (50~100 µg/ml) plate for colony selection.

Recommended colony PCR condition


(SMOBIO's TP1200 ExcelTaq™ 5X PCR Master Dye Mix is suggested)

Template	Single colony
Forward primer	0.1 – 0.5 µM
Reverse primer	0.1 – 0.5 µM
5X PCR Master Dye Mix	5 µl
ddH ₂ O	Variable
<hr/>	
Total volume	25 µl

94°C	2 min	} 30~35 cycles
94°C	30 sec	
50°C	30 sec	
72°C	30 sec/kb	
72°C	1 min	

Note: The amplicon size of the colony PCR is “insert size +152 bp” when using pGet-For Primer and pGet-Rev Primer.

The plasmid map and cloning sites of pGet1.1 Vector



pGet-For primer

5'-GGG ATC CTC GAA GTT AAA GAT GAT TAC GGT GAA TTC AGA ATT CTA CTA GTG GCA GAA
 3'-CCC TAG GAG CTT CAA TTT CTA CTA ATG CCA CTT AAG TCT TAA GAT GAT CAC CGT CTT

GCT AAA CAT CAG GGA AAG GAT ³³⁴ **Cloning site** ³⁹⁵ ATC AAT ATA CGT AAT AAT GGG ATC CTA GTG
 CGA TTT GTA GTC CCT TTC CTA TAG TTA CAT GCA TTA TTA CCA CCC TAG GAT CAC

GGT AAG AGA GGA GAC CAA GAT TTG ATG GCT GCA GGA AAT GCT ATC GAG AGA-3'
 CCA TTC TCT CCT CTG GTT CTA AAC TAC CGA CGT CCT TTA CGA TAG CTC TCT-5'

pGet-Rev primer

Genetic elements of pGet1.1 Vector

Element	Function	Position (bp)
Lethal gene	For screening against self-ligation	1...834
MCS	Multiple cloning site	278...412
Insertion site	The ligation site of blunt end insert	354...355
M13 Rev	M13 reverse primer	838...854
pUC-ori	Initiation of replication	1022...1610
P	The promoter for expressing the ampicillin resistance and lethal gene	2011...2115
Amp ^R	Ampicillin resistance gene	2115...2973
M13 For	M13 forward primer	2975...2990
Primer position		
pGet-For primer	Sequencing of insert, colony PCR	284-305
pGet-Rev primer	Sequencing of insert, colony PCR	415-435

The restriction enzyme with one or two restriction sites on pGet1.1

Enzyme	Positions	Enzyme	Positions
AclI	2302, 2675	FspI	2681
ApaI	1353, 2232	HindIII	204
Avall	2540, 2762	HpaI	568, 822
BamHI	279, 375	NspI	4, 1671
BanI	1923, 2951	PciI	1667, 2995
BglII	434	PssI	1935
Bme1580I	1357, 2236	PstI	418
BsaAI	368, 791	PvuI	2535
BsaHI	2364	PvuII	1847
BseYI	1363	Scal	2423
BspHI	947, 2064	SnaBI	368
BsrDI	2670, 2844(c)	SpeI	322
BsrFI	2819	SspI	2099
CdiI	2351, 2652 (c)	StyI	718
EcoRI	315	TatI	2421

For more details of the sequence of the GetClone™ vector, visit our website www.smobio.com

Related Products

CK1000	Champion E. coli Transformation Kit
CV1100	GetClone PCR Cloning Vector II, 20 RXN
DM2300	ExcelBand 100 bp+3K DNA Ladder, 500 μ l
DM3200	ExcelBand 1 KB Plus (0.1-10 kb) DNA Ladder, 500 μ l
NS1000	FluoroVue Nucleic Acid Gel Stain (10,000 \times), 500 μ l
PM2500	ExcelBand 3-color Regular Range Protein Marker, 250 μ l \times 2
PM2800	ExcelBand 3-color Extra Range Protein Marker, 250 μ l \times 2
RP1000	ExcelRT Reverse Transcriptase, 20,000 U
TF1000	SMO-HiFi DNA Polymerase, 100 U
TF2000	Q-HiFi DNA Polymerase, 100 U
TF3000	G-HiFi DNA Polymerase, 100 U
TP1000	ExcelTaq DNA Polymerase, 500 U
TP5000	ExcelTaq Hot Start II DNA Polymerase, 500 U
TP1200	ExcelTaq 5 \times PCR Master Dye Mix, 200 RXN
TP1260	ExcelTaq 5 \times Fluorescent PCR Master Mix Kit, 200 RXN
VE0100	B-BOX Blue Light LED epi-illuminator, AC 100-240V, 50/60Hz
WM1000	YesBlot Western Marker I, 250 μ l