



SMOBIO

Small Bio, Smart Tool

www.smobio.com

Product information

FluoroBand 50 bp Fluorescent DNA Ladder, Ready to Use

DM1160 **54 μ g / 500 μ l**

Size: 500 μ l

5 μ l per well for 100 gel applications

Storage

Protected from light

Room temperature \geq 6 months

4°C \geq 12 months

-20°C \geq 24 months

Description

The DM1160 FluoroBand 50 bp Fluorescent DNA Ladder is a ready-to-use DNA ladder, which is premixed with high sensitivity DNA binding fluorescent dye and loading dye for direct loading on gel. The DM1160 DNA Ladder is composed of 17 individual DNA fragments: 1.5k, 1.2k, 1k, 900, 800, 700, 600, **500**, 450, 400, 350, 300, 250, **200**, 150 and 100 base pairs derived from a mixture of PCR product and specifically digested plasmid DNA; these bands can be visualized when illuminated with 470 nm blue light or UV light. This product contains two enhanced bands (500 and 200 bp) for easy referencing. In addition, Orange G tracking dye which mimics the migration of 50 bp dsDNA during 1% agarose gel electrophoresis are also added for real time monitoring. Real time observation of the electrophoresis is also possible if compatible light source is fitted to the electrophoresis tank.

Buffering


10 mM Tris-HCl (pH 8.0), 10 mM Na-EDTA, 0.025% Orange G and 10% Glycerol.

Quality Control

1. 17 defined bands are formed during the electrophoresis on the agarose gel.
2. No nuclease activity is observed by a direct nuclease activity assay.

The ladder is supplied in gel loading buffer and is ready to use. Do not heat or dilute before loading.

DM1160



2% Agarose Gel
0.5x TAE Buffer

Note: This DNA Ladder was not designed for precise quantification of DNA mass, but can be used for approximating the mass of DNA by intensity comparison.




SMOBio
Small Bio, Smart Tool

FluoroBand 50 bp Fluorescent DNA Ladder

DM1160

DNA Mass
(ng/ 5 μ l)

Base Pairs
(bp)



2% Agarose Gel
0.5x TAE Buffer



Other information

SMOBiO, Inc. claims all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall SMOBiO, Inc. be liable, whether in contract, tort, warranty, or under any statute or any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including but not limited to the use thereof.

Caution: Not intended for human or animal diagnostic or therapeutic uses.

Related Products

- DM1100 ExcelBand 50 bp DNA Ladder, 500 μ l
- DM2100 ExcelBand 100 bp DNA Ladder, 500 μ l
- DM2160 FluoroBand 100 bp Fluorescent DNA Ladder, 500 μ l
- DM2300 ExcelBand 100 bp+3K Fluorescent DNA Ladder, 500 μ l
- DM2360 FluoroBand 100 bp+3K Fluorescent DNA Ladder, 500 μ l
- DM3100 ExcelBand 1 KB (0.25-10 kb) DNA Ladder, 500 μ l
- DM3160 FluoroBand 1 KB (0.25-10 kb) Fluorescent DNA Ladder, 500 μ l
- DM3200 ExcelBand 1 KB Plus (0.1-10 kb) DNA Ladder, 500 μ l
- DM3260 FluoroBand 1 KB Plus (0.1-10 kb) Fluorescent DNA Ladder, 500 μ l
- DM4100 ExcelBand XL 25 kb DNA Ladder, Broad Range (up to 25 kb), 500 μ l
- DM4160 FluoroBand XL 25 kb Fluorescent DNA Ladder, Broad Range (up to 25 kb), 500 μ l
- DL5000 FluoroDye DNA Fluorescent Loading Dye (Green, 6X), 1ml
- DS1000 FluoroStain DNA Fluorescent Staining Dye (Green, 10,000X), 500 μ l
- VE0100 B-BOXTM Blue Light LED epi-illuminator



B-BOX™ Blue Light LED epi-illuminator

For research use only
2013 ver. 1.2.0