



SMOBIO

Small Bio, Smart Tool

www.smobio.com

Product Information

G-HiFi™ DNA Polymerase

TF3000	100 units	
	G-HiFi™ DNA Polymerase	100 µl
	5X G-HiFi™ Buffer	1200 µl
	dNTPs Mix (2 mM each)	600 µl

Storage

-20°C ≥ 24 months

Description

The G-HiFi™ DNA Polymerase is a new genetically modified, recombinant DNA polymerase suitable for GC-rich templates that are difficult to amplify. The fidelity of G-HiFi™ DNA Polymerase is 70 times higher than that of *Taq* DNA polymerase. The high extension rate of G-HiFi™ DNA Polymerase is achieved by blending the DNA polymerase with an elongation enhancer. The optimized 5X G-HiFi™ Buffer includes special ingredients that suppress non-specific amplification as well as plateau effect produced by conventional PCR. With the optimized 5X G-HiFi™ Buffer, G-HiFi™ DNA Polymerase is capable to amplify most templates, such as longer targets (up to 40 kb from lambda DNA) and that contain GC-rich sequences.

Applications

- Long range PCR amplification for next-generation DNA sequencing
- Generates blunt end amplicons for cloning with GetClone™ PCR cloning vector
- Amplification of GC-rich templates

Features

- 5'→3' DNA polymerase activity
- 3'→5' exonuclease (proofreading) activity
- Suitable for GC-rich templates
- High reaction rate: 7 seconds/kb
- High fidelity: 70 times higher than *Taq* polymerase
- Generates blunt end amplicons
- Vast elongation capability (up to 40 kb)
- Thermo-stable for more than 10 hrs at 95°C.

Storage Buffer

50 mM Tris-HCl (pH 8.0), 50 mM KCl, 0.1 mM EDTA, 1 mM DTT, stabilizer, 50% (v/v) glycerol

Unit Definition

One unit is defined as the amount of enzyme that will incorporate 10 nmol of dNTP into acid-insoluble material in 30 minutes at 74°C.

Recommended PCR Condition

Template	1 – 150 ng
Forward primer	0.1 – 0.5 μM *
Reverse primer	0.1 – 0.5 μM *
5X G-HiFi™ Buffer	10 μl
dNTPs (2 mM each)	5 μl
G-HiFi™ DNA Polymerase	1 μl (1 unit)
H ₂ O	to 50 μl
<hr/>	
Total volume	50 μl

*When amplifying products ≥ 10 kb in length, use primers at a final concentration of 0.1 μM each.

Recommended Primer design

For ≤ 10 kb products:

For general amplification, select primers with a T_m value of $\geq 55^\circ\text{C}$. 20- to 25-mer primers are suitable, or greater than 25-mer in length may provide optimal results.

For >10 kb products:

Select primers with a T_m value of $\geq 65^\circ\text{C}$. 25- to 35-mer primers are suitable. Avoid high GC-content at the 3' end of each primer.

Recommended PCR Program

For GC-rich templates:

98°C	2 min	} 25 ~ 40 cycles
98°C	10 sec	
68°C	10-30 sec/kb	

For ≤ 10 kb products:

98°C	2 min	} 25 ~ 40 cycles
98°C	10 sec	
50~68°C*	15 sec	
68°C	10-30 sec/kb	
68°C	1 min	

For ≥ 10 kb products:

98°C	10 sec	} 25 ~ 40 cycles
68°C	10-30 sec/kb	

* Optimal PCR condition varies according to primers' thermodynamic properties.

Quality Control

Functional Testing

G-HiFi™ DNA Polymerase is tested for performance in the polymerase chain reaction (PCR) using 1 unit of enzyme to amplify a 665 bp target from 10 pg of tested plasmid DNA. The resulting PCR product is visualized as a single band on an ethidium bromide-stained agarose gel.

Nuclease Assay

No contaminating endonuclease or exonuclease activity was detected using pUC19 incubated with G-HiFi™ DNA Polymerase for 4 hours at 37°C.

Residual Nucleotides Assay

No contaminating residual nucleotides were detected from purified G-HiFi™ DNA Polymerase by PCR assay.

Other Information

SMOBIO Technology, Inc. claims all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall SMOBIO Technology, Inc. be liable, whether in contract, tort, warranty, or under any statute or any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including but not limited to the use thereof.

Caution: Not intended for human or animal diagnostic or therapeutic uses.

Related Products

CV1000	GetClone PCR Cloning Vector, 20 RXN
CV1100	GetClone PCR Cloning Vector II, 20 RXN
DM2300	ExcelBand 100 bp+3K DNA Ladder, 500 μ l
DM3100	ExcelBand 1 KB (0.25-10 kb) DNA Ladder, 500 μ l
DM4100	ExcelBand XL 25 kb DNA Ladder, Broad Range (up to 25 kb), 500 μ l
DL5000	FluoroDye DNA Fluorescent Loading Dye (Green, 6 \times), 1 ml
DS1000	FluoroStain DNA Fluorescent Staining Dye (Green, 10,000 \times), 500 μ l
NS1000	FluoroVue Nucleic Acid Gel Stain (10,000X), 500 μ l
TP1000	ExcelTaq DNA Polymerase, 5 U/ μ l, 500 U \times 1
TP1200	ExcelTaq 5 \times PCR Master Dye Mix, 200 RXN
TP1260	ExcelTaq 5 \times Fluorescent PCR Master Mix, 200 RXN
TP2000	ExcelTaq Blood Direct DNA Polymerase, 5 U/ μ l, 500 U
TF1000	SMO-HiFi™ DNA Polymerase, 1U/ μ l, 100 U
TF2000	Q-HiFi™ DNA Polymerase, 1U/ μ l, 100 U
TF3100	2X G-HiFi™ PCR Master Mix, 50 RXN
VE0100	B-BOX™ Blue Light LED epi-illuminator, AC 100-240V, 50/60Hz