



SMOBIO

Small Bio, Smart Tool

www.smobio.com

Product Information

**ExcelBand 3-color Pre-Stained Protein Marker,
High Range (9-245 kDa)**

PM2600

Size: 250 μ l \times 2
5 μ l per well for 100 applications

Storage

4°C \geq 3 months
-20°C \geq 24 months

Loading Volume

- 3 μ l or 5 μ l per loading for clear visualization during electrophoresis on 15-well or 10-well mini-gel, respectively.
- 1.5~2.5 μ l per well for general Western transferring.
- Apply more for thicker (> 1.5 mm) or larger gel.

Description

The PM2600 3-color Pre-Stained Protein Marker is a ready to use three-color protein standard with 12 pre-stained proteins covering a wide range of molecular weights from 10 to 245 kDa in Tris-Glycine buffer (9 to 235 kDa in Bis-Tris (MOPS) buffer and 10-235 kDa in Bis-Tris (MES) buffer). Proteins are covalently coupled with different chromophores for easy identification of bands, with two reference proteins carrying enhanced intensity corresponding to a **green** at 25 kDa and **red** at 75 kDa, respectively, as separated on SDS-PAGE (Tris-Glycine buffer).

The PM2600 3-color Pre-Stained Protein Marker is designed for monitoring protein separation during SDS-polyacrylamide gel electrophoresis, verification of Western transfer efficiency on membranes (PVDF, nylon, or nitrocellulose) and for approximating the size of proteins.

The marker is supplied in gel loading buffer and is ready to use. **Do NOT heat, dilute and add reducing agent before loading.**


Contents

Approximately 0.1~0.4 mg/ml of each protein in the buffer (20 mM Tris-phosphate, pH 7.5 at 25°C), 2 % SDS, 0.2 mM Dithiothreitol, 3.6 M Urea, and 15 % (v/v) Glycerol).

Guide for Molecular Weight Estimation (kDa)

Migration patterns of PM2600 3-color Pre-Stained Protein Marker in different electrophoresis conditions are listed below:

Color	TRIS-GLYCINE	BIS-TRIS (MOPS)	BIS-TRIS (MES)
Blue	245	235	235
Blue	180	170	170
Blue	140	130	130
Blue	100	93	93
Red	75	70	72
Blue	60	53	53
Blue	45	41	42
Blue	35	30	30
Green	25	22	23
Blue	20	18	18
Blue	15	14	14
Blue	10	9	10



4-12% BIS-TRIS gel
MES buffer




SMOBIO
Small Bio, Smart Tool

ExcelBand 3-color Pre-Stained Protein Marker

High Range

PM2600



Note. The apparent molecular weight (kDa) of each protein has been determined by calibration against an unstained protein standard; supplemental data should be considered for more accurate adjustment in different electrophoresis conditions.

Quality Control

Under suggested conditions, PM2600 3-color Pre-Stained Protein Marker resolves 12 major bands in 15% SDS-PAGE (Tris-glycine buffer) and after Western blotting to nitrocellulose membrane.

Caution! Not intended for human or animal diagnostic or therapeutic uses.

Related Products

- PM1500 ExcelBand All Blue Regular Range Protein Marker, 250 μ l \times 2
- PM1600 ExcelBand All Blue Regular Range Plus Protein Marker, 250 μ l \times 2
- PM1700 ExcelBand All Blue Broad Range Protein Marker, 250 μ l \times 2
- PM2400 ExcelBand Pink Blue Protein Marker, 250 μ l \times 2
- PM2500 ExcelBand 3-color Regular Range Protein Marker, 250 μ l \times 2
- PM2700 ExcelBand 3-color Broad Range Protein Marker, 250 μ l \times 2
- PM5000 ExcelBand 3-color Pre-stained Protein Ladder, Regular Range, 250 μ l \times 2
- PM5100 ExcelBand 3-color Pre-stained Protein Ladder, High Range, 250 μ l \times 2
- PM5200 ExcelBand 3-color Pre-stained Protein Ladder, Broad Range, 250 μ l \times 2
- PS1000 FluoroStain Protein Fluorescent Staining Dye (Red, 1,000X), 1 ml
- PS1001 FluoroStain Protein Fluorescent Staining Dye (Red, 1,000X), 1 ml \times 5
- DM1160 FluoroBand 50 bp Fluorescent DNA Ladder, 500 μ l
- DM2160 FluoroBand 100 bp Fluorescent DNA Ladder, 500 μ l

- DM2360 FluoroBand 100 bp+3K Fluorescent DNA Ladder, 500 μ l
- DM3160 FluoroBand 1 KB (0.25-10 kb) Fluorescent DNA Ladder, 500 μ l
- DM3260 FluoroBand 1 KB Plus (0.1-10 kb) Fluorescent DNA Ladder, 500 μ l
- DM4160 FluoroBand XL 25 kb Fluorescent DNA Ladder, Broad Range (up to 25 kb), 500 μ l
- DL5000 FluoroDye DNA Fluorescent Loading Dye (Green, 6X), 1ml
- DS1000 ExcelStain DNA Fluorescent Staining Dye (Green, 10,000X), 500 μ l
- DS1001 ExcelStain DNA Fluorescent Staining Dye (Green, 10,000X), 500 μ l x 5
- VE0100 B-BOX™ Blue Light LED epi-illuminator, AC 100-240V, 50/60Hz