



**SMOBio**

Small Bio, Smart Tool

[www.smobio.com](http://www.smobio.com)

## Product Information

ExcelBand™ series

### Enhanced 3-color High Range Protein Marker (9-245 kDa)

#### PM2610

Size: 250  $\mu$ l  $\times$  2 (for 100 applications)

#### PM2611

Size: 250  $\mu$ l  $\times$  10 (for 500 applications)

## Storage

- 4°C for 3 months
- 20°C for 24 months

## Loading Volume

- **3  $\mu$ l** or **5  $\mu$ l** per loading for clear visualization during electrophoresis on a 15-well or 10-well mini-gel.
- **1.5~2.5  $\mu$ l** per well for general Western blotting transfer.
- Apply more for a thicker (> 1.5 mm) or larger gel.

## Description

The PM2610/PM2611 ExcelBand™ Enhanced 3-color High Range Protein Marker is a ready to use three-color protein standard with 12 pre-stained proteins covering a wide range of molecular weights from 10 to 245 kDa in Tris-Glycine buffer (9 to 235 kDa in Bis-Tris (MOPS) buffer and 10-235 kDa in Bis-Tris (MES) buffer). Proteins are covalently coupled with different chromophores for easy identification of bands, with two reference proteins carrying enhanced intensity corresponding to a **green** at 25 kDa and **red** at 75 kDa, respectively, as separated on SDS-PAGE (Tris-Glycine buffer).

The ExcelBand™ Enhanced 3-color High Range Protein Marker is designed for monitoring protein separation during SDS-polyacrylamide gel electrophoresis, verification of Western transfer efficiency on membranes (PVDF, nylon, or nitrocellulose) and for approximating the size of proteins.


The marker is supplied in gel loading buffer and is ready to use. **Do NOT heat, dilute, or add reducing agents before loading.**

## Contents

Approximately 0.2~0.6 mg/ml of each protein in the buffer (20 mM Tris-phosphate (pH 7.5 at 25°C), 2 % SDS, 3.6 M Urea, and 15 % (v/v) Glycerol).

## Guide for Molecular Weight Estimation

Migration patterns and approximate MWs (kDa) of PM2610/PM2611 ExcelBand™ Enhanced 3-color High Range Protein Marker in different electrophoresis conditions are listed below:

Band	Color	TRIS-GLYCINE	BIS-TRIS (MOPS)	BIS-TRIS (MES)	
1	Blue	245	235	235	
2	Blue	180	170	170	
3	Blue	140	130	130	
4	Blue	100	93	93	
5	Red	75	70	72	
6	Blue	60	53	53	
7	Blue	45	41	42	
8	Blue	35	30	30	
9	Green	25	22	23	
10	Blue	20	18	18	
11	Blue	15	14	14	
12	Blue	10	9	10	




**SMOBio**  
Small Bio, Smart Tool

**ExcelBand™**

**Enhanced 3-color High Range Protein Marker**

**PM2610**



**Note.** The apparent molecular weight (kDa) of each protein has been determined by calibration against an unstained protein standard; supplemental data should be considered for more accurate adjustments in different electrophoresis conditions.



## **Quality Control**

Under suggested conditions, PM2610/PM2611 ExcelBand™ Enhanced 3-color High Range Protein Marker resolves 12 major bands in SDS-PAGE (Tris-Glycine, MOPS, and MES buffer) and after Western blotting transfer to a nitrocellulose membrane.

## **Other information**


SMOBIO Technology, Inc. claims all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall SMOBIO Technology, Inc. be liable, whether in contract, tort, warranty, or under any statute or any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including but not limited to the use thereof.

Caution: Not intended for human or animal diagnostic or therapeutic uses.

## Related Products

- PM1500 ExcelBand All Blue Regular Range Protein Marker, 250  $\mu\text{l}$   $\times$  2
- PM1600 ExcelBand All Blue Regular Range Plus Protein Marker, 250  $\mu\text{l}$   $\times$  2
- PM1700 ExcelBand All Blue Broad Range Protein Marker, 250  $\mu\text{l}$   $\times$  2
- PM2400 ExcelBand Pink Blue Protein Marker, 250  $\mu\text{l}$   $\times$  2
- PM2500 ExcelBand 3-color Regular Range Protein Marker, 250  $\mu\text{l}$   $\times$  2
- PM2510 ExcelBand Enhanced 3-color Regular Range Protein Marker, 250  $\mu\text{l}$   $\times$  2
- PM2600 ExcelBand 3-color High Range Protein Marker, 250  $\mu\text{l}$   $\times$  2
- PM2700 ExcelBand 3-color Broad Range Protein Marker, 250  $\mu\text{l}$   $\times$  2
- PM2800 ExcelBand 3-color Extra Range Protein Marker, 250  $\mu\text{l}$   $\times$  2
- PM5000 ExcelBand 3-color Pre-stained Protein Ladder, Regular Range, 250  $\mu\text{l}$   $\times$  2
- PM5100 ExcelBand 3-color Pre-stained Protein Ladder, High Range, 250  $\mu\text{l}$   $\times$  2
- PM5200 ExcelBand 3-color Pre-stained Protein Ladder, Broad Range, 250  $\mu\text{l}$   $\times$  2

- PS1000 FluoroStain Protein Fluorescent Staining Dye (Red, 1,000X), 1 ml
- DM1100 ExcelBand 50 bp DNA Ladder, 500  $\mu$ l
- DM2200 AccuBand 100 bp DNA Ladder II, 500  $\mu$ l
- DM3100 ExcelBand 1 KB (0.25-10 kb) DNA Ladder, 500  $\mu$ l
- DM5100 ExcelBand Super Range DNA Ladder, (50 bp- 25 kb), 500  $\mu$ l
- VE0100 B-BOX™ Blue Light LED epi-illuminator, AC 100-240V, 50/60Hz
- WM1000 YesBlot Western Marker I, 250  $\mu$ l



YesBlot™ Western Marker I