



SMOBIO

Small Bio, Smart Tool

www.smobio.com

Product Information

ExcelRT Reverse Transcriptase

RP1000 **20,000 units**

ExcelRT Reverse Transcriptase 100 μ l

5X RT Buffer 1 ml

Storage

-20°C \geq 24 months

Description

The ExcelRT Reverse Transcriptase is a recombinant Moloney Murine Leukemia Virus (M-MLV) reverse transcriptase – a RNA dependent DNA polymerase capable of generating the first strand cDNA using RNA template. Designed to reduce RNase H activity and better thermal stability, the ExcelRT Reverse Transcriptase can synthesize routinely first strand cDNA >8 kb at 37~50°C.

Features

- High yield
- Thermostable, up to 50°C, during first strand synthesis
- High processivity, generating cDNA up to 8 kb
- Reduced RNase H ribonuclease activity

Application

- Generation of first strand cDNA from total RNA or mRNA.
- Processivity up to 8 kb.

Storage Buffer

20 mM Tris-HCl (pH 7.5), 200 mM NaCl, 0.1 mM EDTA, 1 mM DTT, stabilizer, 50% (v/v) glycerol

5X RT Buffer

250 mM Tris-HCl (pH 8.3 at 25°C), 50 mM DTT, 375 mM KCl and 15 mM MgCl₂

Unit Definition

One unit is defined as the amount of enzyme that will incorporate 1 nM of dTTP into acid-insoluble material in 10 minutes at 37°C using Poly(A)•oligo(dT)₂₅ as template-primer.

First Strand Synthesis Condition

1. Denature (Mixture A):

Total RNA X μ l (1pg~2 μ g)

Primers 100 μ M d(T)₂₀

or 10 μ M specific primers

or 100 ng/ μ l random hexamers

0.5 μ l

dNTP 1.0 mM (each)

(*High concentration of stock dNTP (10 mM, each) is recommended, allowing greater RNA volume be added in the event of low RNA yield.*)

DEPC-treated H₂O to 10 μ l final vol.

Mix well; incubate at 70°C/5 minutes

Place on ice for at least 1 minute

2. First Strand DNA Buffer (Mixture B) per reaction:

(Master mix can be prepared during or before the denaturing step)

5X RT reaction Buffer 4 μ l

DEPC-treated H₂O 4 μ l

RNase inhibitor 1 μ l

ExcelRT reverse Transcriptase 1 μ l

Final volume 10 μ l

First Strand Synthesis Condition (continued)

3. First Strand cDNA synthesis:

Mixture A (RNA + Primers + dNTP) 10 μ l

Mixture B (First Strand DNA Buffer) 10 μ l

Final Volume 20 μ l

Incubate (25°C/10 minutes)*
37~50°C/50 minutes

4. Termination: 85°C/5 minutes

Keep at 4°C

5. RNA removal[#]: add 1 μ l RNase H into each reaction

37°C/20 minutes

Store cDNA at -20°C or for immediate PCR reaction

*For random hexamers, additional 10 minutes incubation at 25°C is suggested.

[#]Optional Step recommended for long range RT-PCR reaction.

Recommended PCR Condition

(SMOBiO's TP1000 ExcelTaq DNA Polymerase)

cDNA	2~10 μ l
Forward primer	0.1 – 0.5 μ M
Reverse primer	0.1 – 0.5 μ M
10 \times <i>Taq</i> buffer	5 μ l
dNTPs	0.2 mM each
<i>Taq</i> enzyme	0.25 μ l (1.25units)
H ₂ O	to 50 μ l
<hr/>	
Total volume	50 μ l

Recommended PCR Program

94°C	2 min	} 25 ~ 40 cycles
94°C	30 sec	
50~68°C**	30 sec	
72°C	30 sec/kb	
72°C	1 min	

**Optimal PCR condition varies according to primers' thermodynamic properties.

Related Products

DM1100	ExcelBand 50 bp DNA Ladder, 500 μ l
DM2100	ExcelBand 100 bp DNA Ladder, 500 μ l
DM2300	ExcelBand 100 bp+3K DNA Ladder, 500 μ l
DM3100	ExcelBand 1 KB (0.25-10 kb) DNA Ladder, 500 μ l
DM3200	ExcelBand 1 KB Plus (0.1-10 kb) DNA Ladder, 500 μ l
DM4100	ExcelBand XL 25 kb DNA Ladder, Broad Range (up to 25 kb), 500 μ l
DL5000	FluoroDye DNA Fluorescent Loading Dye (Green, 6 \times), 1 ml
DS1000	FluoroStain DNA Fluorescent Staining Dye (Green, 10,000 \times), 500 μ l
TF1000	ExcelTaq SMO-HiFi DNA Polymerase, 1 U/ μ l, 100 U \times 1
TP1000	ExcelTaq DNA Polymerase, 500 U \times 1
TP1200	ExcelTaq 5 \times PCR Master Dye Mix, 200 RXN
TP1260	ExcelTaq 5 \times Fluorescent PCR Master Mix Kit, 200 RXN
TP2000	ExcelTaq Blood Direct DNA Polymerase, 5 U/ μ l, 500 U \times 1
TP2100	ExcelTaq Blood Direct PCR Master Mix Kit, 200 RXN



SMOBio
Small Bio, Smart Tool