



SMOBIO

Small Bio, Smart Tool

www.smobio.com

Product Information

ExcelRT™ series

One-Step RT-PCR Kit

RP1100 50 RXN

Taq /RT Enzyme Mix

50 µl

2X One-step buffer

750 µl x 2

Storage

-20°C for 24 months

Features

- Generates dsDNA directly from RNA samples
- High sensitivity and yields
- High reproducibility, reduced pipetting errors
- High throughput
- Reduced RNase H ribonuclease activity

Description

The ExcelRT™ One-Step RT-PCR Kit is designed for the reverse transcription and PCR amplification of a specific target RNA from either total RNA or mRNA. The ExcelRT™ One-Step RT-PCR Kit provides the user an alternative to the lengthy two step processes (first strand generation and amplification) by using a single mixture, single tube, one step reaction. The ExcelRT™ One-Step RT-PCR Kit contains a 2X reaction premix consisting of an optimized buffer, dNTPs, Mg²⁺ and enzyme stabilizer, and a blend of recombinant reverse transcriptase and *Taq* DNA polymerase. The ExcelRT™ One-Step RT-PCR Kit allows the user to complete the RT-PCR process using a thermocycler in a single reaction setting. The ExcelRT™ One-Step RT-PCR Kit is capable of detecting even trace amounts of target RNA and ideal for target RNA amplification and analysis

General Procedure:

1) Prepare the reaction mix by combining the indicated volume of components of One-Step RT-PCR:

2)

2X One-Step Buffer	25 μ l
Upstream primer (10 μ M) *	1 μ l
Downstream primer (10 μ M) *	1 μ l
Taq /RT enzyme mix	1 μ l
RNA template **	X μ l (10 pg~1 μ g)
DEPC-Treated H ₂ O	Up to 50 μ l
<hr/>	
Final volume	50 μ l

3) Place reaction tubes directly into thermocycler and initiate first strand synthesis – amplification reaction.

* RP1100 is optimized for sequence specific primers.

** Adding RNA template last is recommended.

4) Thermal cycler program

<u>Reaction Steps</u>	<u>Temperature and Time</u>	
Preheat	until block temperature reaches 42~50°C	
<u>cDNA synthesis</u>		
Reverse transcription	42~55°C	30 min
Denature	94°C	2 min
<u>Amplification</u>		
Denature	94°C	20 sec
Annealing [#]	50~68°C	30 sec
Extension	72°C	30 sec/kb
} 30~35 Cycles		
<u>Final extension</u>		
Extension	72°C	1 min

[#]Optimal PCR condition varies according to primers' thermodynamic properties.

Note: RP1100 is not designed for multiplexing DNA with internal control, due to the competition of highly abundant housekeeping gene mRNA.

Other Information

SMOBIO Technology, Inc. claims all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall SMOBIO Technology, Inc. be liable, whether in contract, tort, warranty, or under any statute or any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including but not limited to the use thereof.

Caution: Not intended for human or animal diagnostic or therapeutic uses.

Related Products

CK1000	Champion E. coli Transformation Kit
RI1000	RNAok RNase Inhibitor, 2000 U
RP1000	ExcelRT Reverse Transcriptase, 20,000 U
RP1300	ExcelRT Reverse Transcription Kit, 100 RXN
RP1400	ExcelRT Reverse Transcription Kit II, 100 RXN
DM1100	ExcelBand 50 bp DNA Ladder, 500 μ l
DM2300	ExcelBand 100 bp+3K DNA Ladder, 500 μ l
DM3100	ExcelBand 1 KB (0.25-10 kb) DNA Ladder, 500 μ l
DM4100	ExcelBand XL 25 kb DNA Ladder, Broad Range (up to 25 kb), 500 μ l
DL5000	FluoroDye DNA Fluorescent Loading Dye (Green, 6 \times), 1 ml
NS1000	FluoroVue Nucleic Acid Gel Stain (10,000X), 500 μ l
TQ1100	ExcelTaq 2X Q-PCR Master Mix (SYBR, no ROX), 200 RXN
TQ1110	ExcelTaq 2X Q-PCR Master Mix (SYBR, ROX), 200 RXN
TQ2110	ExcelTaq 2X Q-PCR Master Mix (TaqMan, ROX), 200 RXN

Related Products, continued

TF1000	SMO-HiFi DNA Polymerase, 100 U
TF2000	Q-HiFi DNA Polymerase, 100 U
TF3000	G-HiFi DNA Polymerase, 100 U
TP1000	ExcelTaq Taq DNA Polymerase, 500 U
TP1200	ExcelTaq 5× PCR Master Dye Mix, 200 RXN
TP5000	ExcelTaq Hot Start II DNA Polymerase, 500 U
TP2100	ExcelTaq Blood Direct PCR Master Mix Kit, 200 RXN
VE0100	B-BOX Blue Light LED epi-illuminator, AC 100-240V, 50/60Hz



B-BOX™ Blue Light LED epi-illuminator



SMOBio
Small Bio, Smart Tool

For research use only
2017 ver. 2.1.0